

Simplifying airway care.

Single-use respiratory devices for intensive care, operating rooms and anaesthesia.


• INTENSIVE CARE

Made for the ventilated patient.

One patient, one device. Easy to handle at the bedside, shift after shift.



Single-use · ICU and anaesthesia



1 in 6

ventilated ICU patients in Europe
develop pneumonia on the
ventilator.




2x
odds of death

€21k
extra cost per case

+8.5 days
on the ventilator

OUR MISSION

Make life easier. For nurses and patients.



B-Care
Oral Care

24h Oral Hygiene Kit



Cuffix

Cuff Pressure Regulator




24H ORAL HYGIENE KIT

B-Care Oral Care

● TECHNICAL DATA

Brush	360 degree, soft TPE
Fluid	4 ml pipette
Dosing	Thumb control
Suction	Integrated
Set	Toothbrush, 2 pipettes (4 ml), swab, saliva ejector
Use	Single use, adults

**360° soft brush**

360° soft, rubber-like (TPE) brush: easy access to the patient's mouth

**Fluid in the handle**

Disinfectant-moisturizing fluid stored in handle, dispensed via brush tip

**24-hour protocol**

Simple tracking of 24-hour oral hygiene protocol

**One hand**

Single-hand fluid usage

[View details](#)

CUFF PRESSURE REGULATOR

Cuffix

All-In-One Pressure Regulation

● TECHNICAL DATA

Range	20 to 30 cmH ₂ O
Runtime	Up to 10 days
Power	Internal 3 V battery
Control	Passive
Line	1.5 m extension
Use	Single patient



Reads and regulates

Continuous, precise cuff pressure reading & automated regulation



Stable pressure

Preventing ETT cuff pressure fluctuations and inflation risks



Visible

Light indicators for distance assessment



Fast

Very quick pressure equalization



Easy

Quick and easy handling



Mobile

Disposable and mobile device

View details →

WHAT WE ARE WORKING ON


Closed Suction System.

Suction and wiping of the tube's inner wall with a balloon, assisted by fluid.

All with one push of a button, without changing the care routine.

Coming soon

CLOSED SUCTION SYSTEM · PREVIEW



THE COMPANY

We are AWM.

- **ISO 13485**

Certified quality management. 2026 audit without findings.

- **Warsaw Stock Exchange**

Listed since 2014, ticker AWM.

- **CE under MDR**

CE marked under the EU Medical Device Regulation.

- **Patented**

B-Care, Cuffix and CSS are protected by patents.



Headquarters · ul. Szara 10, 00-420 Warsaw, Poland

One dedicated partner per market.

We work with one distributor in each market. You talk to us, you sell, we support you.



1

Talk to us.

Tell us about your market and the hospitals you serve.

2

Agree on the terms.

Volumes and conditions for your market, set together. One partner per market: in yours, that is you.

3

Get your support material.

Registration documents, IFU, tutorials, trailer and one-pagers, ready to use.

Let's simplify care in your country.

We work with one distributor in each market. You talk to us, you sell, we support you.

sales@airway-medix.com


airway-medix.com

REVOLUTIONIZING MEDICINE, SIMPLIFYING CARE




B-CARE ORAL CARE
AIRWAY MEDIX S.A.

EXPLODED VIEW
2 PARTS + PIPETTE · NOT TO SCALE



[◀ Back to B-Care](#)



CUFFIX
AIRWAY MEDIX S.A.

EXPLODED VIEW
6 PARTS · NOT TO SCALE

[◀ Back to Cuffix](#)